





P	Code	Description	Q	Assembly prescription	Validity	P	Code	Description	Q	Assembly prescription	Validity
1	1560580	Plug	1								
2	540290	Screw	TCEI M8x25	1	C=25Nm						
3	1660210	Nut	M8	1							
4	1560400	Knob		1							
5	1080070	Pin	stainless steel	1							
6	1560420	Screw		1	Lubricate with grease Molykote PG54						
7	2260330	Guide		1	C=35Nm / Loxeal 55-14						
8	2260140	Wobble plate		1							
9	1560350	Spring		1							
10	2260340	Wobble plate		1							
11	2260060	Piston		1							
12	2260100	O-ring	Ø 6,02x2,62	2	Lubricate with grease Molykote PG54						
13	2260040	Piston guide		1							
14	740290	O-ring	Ø 14x1,78	1	Lubricate with grease Molykote PG54						
15	394280	O-ring	Ø 12,42x1,78	1	Lubricate with grease Molykote PG54						
16	2761290	Piston	lower	1	C=3Nm / Loxeal 83-21						
17	2260080	Seat		1							
18	1260740	O-ring	Ø 10,77x2,62	1	Lubricate with grease Molykote PG54						
19	2260010	Valve body		1							
20	1560151	O-ring	Ø 5,94x3,53	1	90 Sh						
21	2260030	Shutter		1							
22	1560140	Spring		1							
23	770130	O-ring	Ø 20,35x1,78	1							
24	2260020	Fitting	3/8" G	1	C=40Nm						
F	2876	Gasket kit		1							

Revision: 1

Oil type: -

Validity: < XX/YY = From week/year

Date: 22/04/2025

Oil quantity: -

C= Tightening torque (Tolerance +0 ÷ -10%)

* = Not present for electric pump